

Malaria Detection

DEEP LEARNING PROJECT

MIT Applied Data Science Program

Capstone Project: Deep learning

Delia McNamara

14 August 2024

CONTENTS / AGENDA

- Executive Summary
- Problem Definition
- Solution summary
- Data Exploration
- Model Performance Comparison
- Final Solution Design Proposal
- Recommendations for implementation

EXECUTIVE SUMMARY

- Use deep learning to diagnose Malaria infection.
 - ◇ Create convolutional neural network (CNN) model to classify imaged cells as parasitized or uninfected.
 - ◇ Intention to utilize the CNN to overcome global deficit of healthcare access in diagnosing Malaria.
- 5 CNNs created and compared.
 - ◇ Focus on overall accuracy & recall for parasitized samples.
 - ◇ 3 CNN models with hidden layers including convolutional, dense, normalization, dropout, and flatten layers.
 - ◇ 1 CNN model with data augmentation.
 - ◇ 1 pre-trained CNN model (VGG16).
- Best model (CNN Model 3) had an accuracy of 97% and a recall for parasitized samples of 98%.
 - ◇ Tried this model with data augmentation.
 - Worse performance with augmentation.
- Other models ranged from an accuracy of 93.1 to 98.0%
- CNN Model 3 recommended for use in diagnostics for overall accuracy, best fit, and lowest proportion of false negatives.

PROBLEM DEFINITION

- 400,00 deaths worldwide caused by Malaria in 2019.
- Malaria is caused the *Plasmodium* parasite, which can be imaged in cells for diagnosis.
- 229 million cases of Malaria diagnosed in 2019.
 - ◇ Traditional diagnosis requires a lab technician to detect parasite presence in sample via microscope.
 - This labor issue limits global access to proper healthcare options for Malaria diagnosis.
 - ◇ Strong need for quick & accurate diagnosis.
- Objective to build CNN models to predict the presence of *Plasmodium* parasite in imaged tissue samples.
 - ◇ Decrease labor needs and increase access to diagnostic tools without compromising patient outcomes.
 - ◇ Focus on high accuracy while minimizing false negatives in models.
 - False negatives are patients who have Malaria, but test results indicate they do not have Malaria.
 - Extremely risky for patients – if they test negative, they will likely not be given the treatment they need.
 - Potentially lethal outcomes for these patients.


SOLUTION SUMMARY

- Data preprocessed at import for consistency of dimensions.
 - ◇ Imported pre-split into test & train sets.
- Data exploration performed to visualize images & how augmentation may affect them.
 - ◇ Images displayed in original RGB format.
 - ◇ Average image calculated & visualized.
 - ◇ Gaussian blur applied.
 - ◇ Images converted to HSV format.
- Labels one-hot encoded for use in neural network.
- CNN Model 1 performed with convolutional, dense, dropout, max pooling, & flatten layers.
- CNN Model 2 performed with similar layers but used Tanh activation rather than ReLU.
- CNN Model 3 performed with similar layers but used Leaky ReLU activation, added normalization, & removed one dropout layer for better fit.
- CNN Model 3 then performed with data augmentation.
 - ◇ Used Image Data Generator to augment.
 - ◇ Horizontal flip, rotation, & zoom applied.
- Pre-trained VGG16 model performed.

DATA EXPLORATION

DATA OVERVIEW

- 27,558 images total in dataset.
- Data labeled, split into train & test sets prior to import.
 - ◇ 24,958 train images.
 - 12,376 images labeled uninfected.
 - 12,582 images labeled parasitized.
 - ◇ 2,600 test images.
 - 1,300 images labeled uninfected.
 - 1,300 images labeled parasitized.



→ Images resized to 64 x 64 px at import for uniformity.

◇ Prevent model from training on imaging artifacts.

→ Cell images fairly uniform.

◇ All have black background.

◇ Parasites appear to have a pink to purple hue.

→ Variable hue for cells in images.


◇ Shades of pink, orange, purple, blue.

→ Some uninfected images appear to have a brightness issue with the center of the cell more brightly lit than the perimeter.


CNN PREPROCESSING

NO IMAGE AUGMENTATION

TRAIN SET



TEST SET



→ Dataset pre-split for model training & validation.


- ◇ 24,958 images to train (90%).
- ◇ 2,600 images to test (10%).
- ◇ 0 represents uninfected, 1 represents parasitized.

→ Good balance between uninfected & parasitized samples.

- ◇ Model will not be biased by the label split.

→ Train & test set values normalized & converted to float data type.

- ◇ Imported with range of pixel values from 0 to 255.
- ◇ Normalized to values of 0 to 1.
- ◇ Facilitates network learning by ensuring pixel values are within a similar range.




→ Calculated mean image for Parasitized & Uninfected.

→ Nearly indistinguishable.


◇ Good for model to learn to distinguish between the two based on the appearance of the parasite and not imaging artifacts or bias.

- E.g. all uninfected samples imaged on a white background, all infected samples imaged on a black background.

WITH IMAGE AUGMENTATION



- Converted train & test images from RGB to HSV.
- Visualize options for data augmentation.
- Images more saturated with less discernable value ranges.
- Consistent cell color between images.
 - ◇ Better than original RGB.
- Parasites appear to have a yellow hue.
- One image contains many blue pixels.
 - ◇ Some anomalies in this image format.



→ Gaussian blur applied to further visualize options for data augmentation.

→ Images less sharp.


→ Parasites still visibly pink to purple.

◇ Edges are less distinguishable.

→ Cell colors still variable between shades of pink, orange, blue, purple.

→ Likely not a good use of augmentation.

◇ Blurring is better suited for reducing features in overly complex images – these images are simple.



→ Image Data Generator used to augment images.

→ Horizontal flip, rotation, & zoom applied.

→ Add variability for the model to learn.


LABEL PREPROCESSING

→ One-hot encoded train & test set labels.

- ◇ 2 labels for each.
- ◇ Allows a dataset originally labeled with categorical labels to be used in neural network.

MODEL PERFORMANCE COMPARISON

CNN MODEL 1



→ CNN model with:

- ◇ 3 convolutional layers.
 - Output size same as input.
 - Kernel size 3 x 3.
 - ReLU activation
- ◇ 3 max-pooling layers.
 - Pool size 2 x 2.
- ◇ 4 Dropout layers
- ◇ 1 flatten layer.
- ◇ 2 dense layers.
- ◇ 1 output layer.
 - 2 nodes.
 - Softmax activation.

→ Compiled with:


- ◇ Loss: categorical cross-entropy.
- ◇ Optimizer: Adam.
 - Learning rate: 0.001.

→ Fitted with:

- ◇ Callback for early stopping with validation loss monitor to prevent further epochs if loss is no longer decreasing.
- ◇ Validation split: 20%.
- ◇ Batch size: 32.
- ◇ Called with 20 epochs.
 - Ran 8 epochs.

→ Model performance on test data:

- ◇ Test accuracy 98.0%



	precision	recall	f1-score	support
0	0.97	0.99	0.98	1300
1	0.99	0.97	0.98	1300
accuracy			0.98	2600
macro avg	0.98	0.98	0.98	2600
weighted avg	0.98	0.98	0.98	2600

→ 33 false negatives, 18 false positives.


- ◇ More images classified as uninfected when actually parasitized.
- ◇ Seek to minimize this issue – false negatives are weighed heavily when the patient outcome can be lethal.

→ Good precision at 97% for uninfected & 99% for parasitized.

- ◇ Few false positives.

→ Good recall at 99% for uninfected, OK at 97% for parasitized.

- ◇ More important metric to avoid lethality to patients.
- ◇ Good result, but can be improved.




→ Validation accuracy higher overall than test accuracy.

◇ May indicate some underfitting.

→ Validation loss stabilizes with training loss by epoch 5.

◇ Underfitting is not too bad.

CNN MODEL 2



→ CNN model with:

- ◇ 3 convolutional layers.
 - Output size same as input.
 - Kernel size 3 x 3.
 - 2 ReLU activation.
 - 1 tanh activation.
- ◇ 3 max-pooling layers.
 - Pool size 2 x 2.
- ◇ 4 Dropout layers
- ◇ 1 flatten layer.
- ◇ 2 dense layers.
 - Tanh activation.
- ◇ 1 output layer.
 - 2 nodes.
 - Softmax activation.

→ Compiled with:


- ◇ Loss: categorical cross-entropy.
- ◇ Optimizer: Adam.
 - Learning rate: 0.001.

→ Fitted with:

- ◇ Callback for early stopping with validation loss monitor to prevent further epochs if loss is no longer decreasing.
- ◇ Validation split: 20%.
- ◇ Batch size: 32.
- ◇ Called with 20 epochs.
 - Ran 6 epochs.

→ Model performance on test data:

- ◇ Test accuracy 96.9%



	precision	recall	f1-score	support
0	0.96	0.98	0.97	1300
1	0.98	0.95	0.97	1300
accuracy			0.97	2600
macro avg	0.97	0.97	0.97	2600
weighted avg	0.97	0.97	0.97	2600

→ 59 false negatives, 22 false positives.

- ◇ More images classified as uninfected when actually parasitized.
- ◇ Seek to minimize this issue – false negatives are weighed heavily when the patient outcome can be lethal.

→ OK precision at 96% for uninfected & 98% for parasitized.



- ◇ Few false positives.

→ Fair recall at 98% for uninfected, poor recall at 95% for parasitized.

- ◇ More important metric to avoid lethality to patients.
- ◇ Poor result for parasitized samples.

→ CNN Model 1 fewer false negatives.

- ◇ Better option for patient outcomes.




→ **Validation accuracy higher overall than training accuracy.**

◇ **May indicate some underfitting.**

→ **Validation loss stabilizes with training loss by epoch 3.**

◇ **Underfitting is not too bad.**

CNN MODEL 3



→ CNN model with:

- ◇ 3 convolutional layers.
 - Output size same as input.
 - Kernel size 3 x 3.
 - Leaky ReLU activation.
- ◇ 3 max-pooling layers.
 - Pool size 2 x 2.
- ◇ 3 Dropout layers
- ◇ 1 flatten layer.
- ◇ 2 dense layers.
 - Leaky ReLU activation.
- ◇ 1 output layer.
 - 2 nodes.
 - Softmax activation.

→ Compiled with:


- ◇ Loss: categorical cross-entropy.
- ◇ Optimizer: Adam.
 - Learning rate: 0.001.

→ Fitted with:

- ◇ Callback for early stopping with validation loss monitor to prevent further epochs if loss is no longer decreasing.
- ◇ Validation split: 20%.
- ◇ Batch size: 32.
- ◇ Called with 20 epochs.
 - Ran 3 epochs.

→ Model performance on test data:

- ◇ Test accuracy 97.0%



	precision	recall	f1-score	support
0	0.98	0.96	0.97	1300
1	0.96	0.98	0.97	1300
accuracy			0.97	2600
macro avg	0.97	0.97	0.97	2600
weighted avg	0.97	0.97	0.97	2600

→ 28 false negatives, 51 false positives.

- ◇ Fewer images classified as uninfected when actually parasitized.

→ Fair precision at 98% for uninfected & 96% for parasitized.


- ◇ Few false negatives, more false positives than previous models.

→ Poor recall at 96% for uninfected, fair recall at 98% for parasitized.

- ◇ More important metric to avoid lethality to patients.
- ◇ Better result for parasitized samples.

→ This model fewest false negatives so far.

- ◇ Better option for patient outcomes.



→ Validation & training accuracy values are similar after 1 epoch.

◇ Good model fit very quickly.

→ Validation loss higher than training loss after 1 epoch.

◇ Some overfitting may be occurring.


→ Fewer epochs.

◇ Favorable for computation, but may be preventing loss from decreasing.

→ Model accuracy is highest for this model out of all performed.

CNN MODEL 3 WITH DATA AUGMENTATION

→ Data augmentation with Image Data Generator noted previously.



→ CNN model with:

- ◇ 3 convolutional layers.
 - Output size same as input.
 - Kernel size 3 x 3.
 - Leaky ReLU activation.
- ◇ 3 max-pooling layers.
 - Pool size 2 x 2.
- ◇ 3 Dropout layers
- ◇ 1 flatten layer.
- ◇ 2 dense layers.
 - Leaky ReLU activation.
- ◇ 1 output layer.
 - 2 nodes.
 - Softmax activation.

→ Compiled with:


- ◇ Loss: categorical cross-entropy.
- ◇ Optimizer: Adam.
 - Learning rate: 0.001.

→ Fitted with:

- ◇ Callback for early stopping with validation loss monitor to prevent further epochs if loss is no longer decreasing.
- ◇ Validation split: 20%.
- ◇ Batch size: 32.
- ◇ Called with 20 epochs.
 - Ran 7 epochs.

→ Model performance on test data:

- ◇ Test accuracy 94.5%



	precision	recall	f1-score	support
0	0.90	1.00	0.95	1300
1	1.00	0.89	0.94	1300
accuracy			0.95	2600
macro avg	0.95	0.95	0.95	2600
weighted avg	0.95	0.95	0.95	2600

→ 138 false negatives, 4 false positives.

- ◇ More images classified as uninfected when actually parasitized.
- ◇ Seek to minimize this issue – false negatives are weighed heavily when the patient outcome can be lethal.

→ Very poor precision at 90% for uninfected, great at 100% for parasitized.

- ◇ Very few false positives for uninfected.

→ Very good recall at 100% for uninfected, very poor recall at 89% for parasitized.


- ◇ More important metric to avoid lethality to patients.
- ◇ Poor result for parasitized samples.

→ Poor F1 scores – precision & recall are not well-balanced.

→ Too many false negatives to use this model.

→ CNN Model 3 fewer false negatives.

- ◇ Better option for patient outcomes.




→ Training & validation accuracies are similar after 2 epochs, diverge slightly, then reconverge after 7 epochs.

◇ Little overfitting.

→ Training & validation loss converge quickly by epoch 2.

→ The fit is good, but accuracy is not very high.

PRE-TRAINED MODEL VGG16



→ CNN model with:

→ Pre-defined layers:

- ◇ 13 convolutional layers.
- ◇ 5 max-pooling layers.

→ Added classification layers:

- ◇ 1 Dropout layer.
- ◇ 1 flatten layer.
- ◇ 2 dense layers.
 - ReLU activation.
- ◇ 1 normalization layer.
- ◇ 1 output layer.
 - 2 nodes.
 - Softmax activation.

→ Compiled with:


- ◇ Loss: categorical cross-entropy.
- ◇ Optimizer: Adam.
 - Learning rate: 0.001.

→ Fitted with:

- ◇ Callback for early stopping with validation loss monitor to prevent further epochs if loss is no longer decreasing.
- ◇ Validation split: 20%.
- ◇ Batch size: 32.
- ◇ Called with 10 epochs.
 - Ran 5 epochs.

→ Model performance on test data:

- ◇ Test accuracy 93.1%



	precision	recall	f1-score	support
0	0.96	0.89	0.93	1300
1	0.90	0.97	0.93	1300
accuracy			0.93	2600
macro avg	0.93	0.93	0.93	2600
weighted avg	0.93	0.93	0.93	2600

→ 43 false negatives, 137 false positives.

◇ More images classified as parasitized when actually uninfected.

→ Poor precision at 96% for uninfected, very poor at 90% for parasitized.

◇ Many false positives for uninfected.

→ Very poor recall at 89% for uninfected, fair recall at 97% for parasitized.



◇ More important metric to avoid lethality to patients.

→ Poor F1 scores – precision & recall are not well-balanced.

→ Fewer false negatives in this model compared to others.

→ CNN Model 3 fewer false negatives.

◇ Better option for patient outcomes.



→ Training & validation accuracies diverge greatly at epoch 1, become more similar at epoch 2 but training remains higher.

◇ Indicates some overfitting.

→ Training & validation loss diverge greatly at epoch 1, become more similar at epoch 2 but validation remains higher.

◇ Indicates overfitting.

→ Overall, one of the poorer fitting models in comparison to the others.

FINAL SOLUTION DESIGN PROPOSAL

→ CNN Model 3 with no data augmentation displayed the best performance for key metrics.

- ◇ Fewest false negatives.

- Extremely important to avoid lethal outcome for patients due to misclassification.

- ◇ Few false positives, but more than false negatives.

- Ideal to minimize, but not at the expense of minimizing false negatives.

- ◇ Great model accuracy at 97%

- Overall predictive power is good.

- ◇ Validation loss indicated some potential overfitting, but accuracy convergence between training & validation sets was contradictory.

RECOMMENDATIONS FOR IMPLEMENTATION

→ Recommend using CNN Model 3 as a tool for Malaria diagnosis.

→ CNN Model 3 displayed the best balance of fit and accuracy.

◇ CNN Model 1 had a higher overall accuracy at 98%.

- Underfitting in training accuracy for this model indicates it may not be able to extrapolate to new images well.

- Fit is important for resiliency in the model to be able to handle images that vary in quality – an issue likely to be seen in areas with limited access to diagnostic resources.

◇ Accuracy of 97% was still one of the highest of all models.

→ CNN Model 3 displayed the lowest proportion of false negatives

◇ One of the key metrics being used to evaluate these models.

- Must avoid diagnosing patients as uninfected if they are parasitized to ensure they are treated for the parasite.